



CLINiC

Quick Start Guide (for Cisco)

Scan for User Guide



Do not manually adjust camera.



**Designed to be powered on 24/7/365
Please do not turn off the unit.**


Powering on/off

Connect the peripherals to the HDMI port(s) of the CLINiC before connecting the CLINiC to AC power.

The system goes to sleep after a period with no activity. It will wake up automatically from an incoming call or pressing any button.

Controls

1. Answer incoming call.
2. End or reject call.
3. Enable video privacy from the CLINiC camera during a call, or revert to normal video transmission.
4. Display and transmit images and data from Horus Scope.
5. Display and transmit images and data from the HDMI/PC port.
6. Zoom camera out.
7. Zoom camera in.
8. Pan and tilt the camera
9. Enable or disable self-view mode
10. Decrease speaker volume.
11. Increase speaker volume
12. Mute or unmute microphone.
13. Enable or disable headphone mode for private listening.
14. Enable or disable stethoscope mode.
15. "On Air" Indicator – Illuminates when system is in an active call.
16. Headphone audio output.
17. Electronic stethoscope audio input.



Accessing the Codec Web Interface

To access the web interface

1. In your web browser address line, enter the system IP address, for example, <http://10.11.12.13>. To display the system IP address, **press and hold the Hang Up button on the CLINiC control panel for 5 seconds or more**, and then release. This will display the system information on the CLINiC screen, including the system IP address.
2. Enter the Admin ID as the user name (default is: admin), and enter the Admin Remote Access Password, if one is set.
3. The Admin screen can now be used to navigate through the menus to the required location to change a setting, as outlined in the following pages.
4. It is important that after any changes to the configuration are made these are saved using the SAVE button in the top right hand corner.



- To enable stand-by mode: Set Standby Control to On and set a Standby delay time
- Set Standby Delay: Input the required number of minutes between 1 and 480
- Save Configuration changes: Select Save

Change Conference Settings

- The default setting is set for the CLINiC to auto answer incoming calls
- This can be changed to manual answering using the connect button on the system control panel
- To change conference call settings, navigate to Setup>Configuration>Conference
- To Disable Auto Answering of Incoming Calls: Set AutoAnswer Mode to Off
- To Enable Auto Answering of Incoming Calls: Set AutoAnswer Mode to On
- Save Configuration changes: Select Save

System Auto-Dial

- A single contact can be auto-dialed from the CLINiC using the connect key
- To add an auto-dial contact, navigate to Setup>Configuration>Facility/Service
- Locate Service 5 and select Call Type: Video
- Input a identifying User name and associated URL
- Select Type: Other
- Select: Save

Audio Levels

You can adjust the level of each audio input transmitted from the codec to optimize the audio experience at the remote end.

- **To adjust individual audio level inputs**, navigate to: Setup>Configuration>Audio (Suggested changes and adjustments screen shot is shown on following page)
- **Microphone 1:** Integrated System Microphone (Default value:18) Adjust as required
Set Mode to: ON
Set Dereverberation to: OFF
Mode to: ON
Noise Reduction to: ON
- **Microphone 2:** Electronic Stethoscope Input (Default value:18) Adjust as required
Set Mode to: OFF
Set Dereverberation to: OFF
Mode to: OFF
Noise Reduction to: OFF
- **Microphone 3:** Auxiliary Microphone Input (Default value:18) Adjust as required
Set Mode to: ON
Set Dereverberation to: OFF
Mode to: ON
Noise Reduction to: ON

Following changes to any levels or settings, Select: SAVE