



MedView Quick Start Guide (for Cisco)

Scan for User Guide



Do not manually adjust camera.



Designed to be powered on 24/7/365 Please do not turn off the unit.

The MedView will power on after being plugged in for 2 minutes. The camera will automatically adjust to a privacy setting.



To reset the camera, either hold down the escape key or press the volume up and down buttons simultaneously for 5 seconds on the remote. If no remote is attached, reset by unplugging for 5 seconds.




The **green phone** icon control is used to make a call and accept an incoming call.



The **red phone icon** control is used to end a call. The system will go into standby mode after a call is dropped and the camera will move into the privacy position.

Controls

1. Pan and tilt the camera (2nd function: moves self-view)
2. Zoom camera in
3. Zoom camera out
4. Increase speaker volume
5. Decrease speaker volume
6. Accept new PIP location
7. ESC - Reset camera
8. Answer incoming call/place call
9. End or reject call
10. Mute/Unmute microphone
11. Enable/Disable video privacy
12. Display/Transmit presentation source
13. Enable/Disable self-view mode



Accessing the Admin Settings

Accessing the Codec Web Interface

To access the web interface:


1. In your web browser address line, enter the system IP address, for example, <http://10.11.12.13>.

The system IP address is always displayed in the bottom left corner of the screen when the CLINiC is on.

2. Enter the Admin ID as the user name (default is admin), and enter the Admin Remote Access Password, if one is set.

3. The Admin screen can now be used to navigate through the menus to the required location to change a setting, as outlined in the following pages.

4. It is important that after any changes to the configuration are made these are saved using the SAVE button in the top right hand corner



System Auto-Dial (Optional Tethered Remote Only)



A single contact can be auto-dialed from the MedView using the connect key.

To add an auto-dial contact, navigate to Setup>Configuration>Facility/Service.

1. Locate: **Service 5**
2. Select Call Type: **Video**
3. Input a identifying: **User name**
4. Input associate: **URL**
5. Select Type: **Other**
6. Select: **Save**

