



JEDMED

Power

Indicator

Horus Scope Series 1

Quick Start Guide

Horus Scope Functions and Components

The Horus Scope functions and device components are listed below.



Basic Horus Scope Operations

To Power on the Horus Scope:

Press and hold the power button at the top of the unit for two seconds.

The Horus Scope will start up. Start-up takes about 20 seconds.

Once the start-up procedure is complete the Horus Scope screen will display the live video feed captured by the lens currently attached to the unit.



To switch between live video and image capture:

1. Press the OK button on the unit to capture images.



2. Press the button to return to live view.